

Date App. Received: _____
Date Order: _____
Date given away: _____



www.capitalcityambucs.org

Amtryke Therapeutic Tricycle Request Form

(Must be filled out completely by adult rider or parent/guardian)

Recipient's Name: _____ Age: _____ Date of Request: _____
Mailing Address: _____ Phone #: _____
City/State/Zip: _____ Email: _____
Diagnosis: _____

If Recipient is Under Age 18

Parent/Guardian Name: _____
If different from above
Mailing Address: _____ Phone #: _____
City/State/Zip: _____ Email: _____

Secondary Contact Name: _____ Phone #: _____

Treating Therapist's Name: _____
Phone #: _____ Email: _____

How did you hear about the Amtryke Therapeutic Tricycle? (Check all that apply)

Therapist Website AMBUCS Member Other: _____

Will you need financial assistance to obtain the tricycle? Yes No

If yes, how much can you pay? _____

Note: Amtryke therapeutic tricycles are distributed based on available funds and need. Individual placements of Amtryke therapeutic tricycles are at the discretion of the local chapter or parent organization.

Tell us about the recipient*: _____

*This information will be made public to help obtain funding. Please don't include information you don't want shared.

★ Including a photo of the recipient will help us obtain a sponsor to help you pay for the Amtryke more quickly. Digital images preferred but we also accept professionally printed glossy photos. No photocopies or folded images.
★ By including a photo, you are giving consent for AMBUCS to use the image online and in print to help obtain a funding.

Signature: _____ Date: _____

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

3/2020

Amtryke Therapeutic Tricycle Waiver Form

(Must be filled out completely by adult rider or parent/guardian)

AMBUCS members nationwide are dedicated to creating opportunities for mobility and independence for people with disabilities by providing Amtryke therapeutic tricycles, offering educational scholarships to therapy students and performing various forms of community service.

Purpose: The Amtryke therapeutic tricycle was designed for people with disabilities. It creates a feeling of freedom, builds self-esteem, strengthens muscles and improves motor coordination and range of motion—all while making exercise fun.

Steering: Initially, the rider may have difficulty turning or changing directions. Encourage the rider to go straight ahead, back up and slowly turn around. On many models there are three steering options for the Amtryke. On the front column of the tricycle you will find two holes for the steering pin. The top hole is straight steering, the bottom hole allows a 20-degree turning radius. Leaving the pin out gives the rider free steering.

Safety Cautions

- Fast speeds and sharp turns can cause the Amtryke therapeutic tricycle to tip or turn over.
- Always wear a helmet when riding an Amtryke. Use of other protective gear is highly recommended.
- Adult supervision required if used by younger or developmentally delayed riders.
- Use caution near vehicles, swimming pools and other bodies of water, hills, alleys and sloped driveways.
- Always wear shoes.
- Never allow more than one rider.
- Use of the steering pin is recommended to prevent over-steering or possible tip-overs.

The information contained in this service is not intended nor implied by National AMBUCS™, Inc. to be professional medical advice by National AMBUCS, Inc. Always seek the advice of your physician, therapist or other qualified healthcare provider prior to starting any treatment or with any question you may have regarding a medical condition. Nothing contained in this document is intended by National AMBUCS, Inc. to be for medical diagnosis or treatment by National AMBUCS, Inc. or on behalf of National AMBUCS, Inc.

In no event shall National AMBUCS, Inc. be liable for any direct, indirect, incidental, consequential, special, exemplary, punitive, or any other monetary or other damages, personal injury or property damages, fees, fines, costs, attorney fees, or liabilities of any kind arising out of or relating in any way to this service or use of the Amtryke® therapeutic tricycle, and/or content or information provided herein.

I agree that the rider listed on the front of this document may be photographed. I also agree that the photos and recipient's name may be used in promotional efforts for National AMBUCS, Amtryke or the local AMBUCS chapter. I further grant AMBUCS the ability to use the photos and name for advertising/publicity purposes without additional compensation, except where prohibited by law. If anybody in my party does not want to be photographed under these same terms, I will let the photographer know as soon as possible.

By signing below, I acknowledge that I have read and understood this liability waiver.

Recipient's Name: _____

Adult Recipient Signature: _____

If Recipient is Under Age 18

Legal Guardian Name: _____

Legal Guardian Signature: _____ Date: _____

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722



PHOTO RELEASE FORM

I, _____, the undersigned agree that my child, _____, may be photographed by National AMBUCS™ Inc. and its subsidiary chapters. I also agree that my daughter/son's photo, name and information may be used in promotional efforts for National AMBUCS™, Inc, AmTryke®, LLC and its subsidiary chapters. I further grant National AMBUCS™, Inc and its subsidiary chapters the ability to use photos name and information for advertising/publicity purposes without additional compensation, except where prohibited by law.

DATE: _____

Witness: _____

Please scan and email this application to:
applications@capitalcityambucs.org

OR

mail to:

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

This form, assessment form and tryke selection form must be received by your local chapter before placement is considered.


Amtryke Assessment Form

(Must be filled out completely by therapist)

Recipient's Name: _____
 Age: _____ Weight (lbs.):* _____ Height (Inches):* _____
 Diagnosis(es):* _____

*This private information is only used to help appropriately fit the rider.

RIDER'S MEASUREMENTS				
Arm Measurements (Inches)				Total Length
Left	A to B:		B to C:	
Right	A to B:		B to C:	
Trunk	A to D:			
Leg Measurements (Inches)				Total Length
Left	D to E:		E to F:	
Right	D to E:		E to F:	



A	Center of Shoulder
B	Center of Elbow
C	Center of Digit Crease
D	Center of Hip
E	Center of Knee
F	Bottom of Foot

Arm Length & Leg Length Measurements are critical to correct Amtryke Selection

Sizing Chart is available online:
www.ambucs.org/riders/wish-list/sizing-chart/

Notes on Provided Measurements (if any):

Dominant Hand:	Left	Right
----------------	------	-------

Helmet Sizing	
Size	Measurement (head circumference)
XXS	18.5" to 19.5"
XS	20.5" to 22"
S/M	22" to 23.6"
L/XL	23.6" to 25.75"

Therapist Name: _____ Is this the treating therapist? Yes No
 Credentials: _____
 Phone: _____ Email: _____
 Facility Name: _____
 Street Address: _____ City: _____ State: _____ Zip: _____
 Is this facility an Amtryke Evaluation Site? Yes No Not Sure

Therapist comments concerning recipient or goals:

This request/assessment is directed to:

- Local AMBUCS Chapter Name: Capital City AMBUCS
 National Wish List (AMBUCS Resource Center)

By signing below, you are signifying that in your professional opinion this rider would benefit from an Amtryke. You assume no liability.

Therapist Signature: _____ Date: _____

Ship Amtryke To			
Name/Facility:	_____	Phone:	_____
Street Address:	_____		
City:	_____	State:	_____ Zip: _____

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
 3207 Dancy St. Austin, Tx 78722

Recipient's Name:

Thanks for choosing an Amtryke therapeutic tricycle!

In order to accommodate the widest variety of people and disabilities, Amtryke offers many tryke models and each can be customized in a variety of ways. Following the steps below will help you choose to the perfect tryke for your client from what might seem like a dizzying array of options.

Remember you can always refer to our website, www.amtrykestore.org, or the Amtryke catalogue for more information and product images.

Step 1: Fill out the first page of the Amtryke Assessment Form.

Step 2: Choose the way the tryke will be propelled: **Hand & Foot, Foot, or Hand**. Your choice should be based on the rider's ability and therapy goals.

Hand & Foot trykes improve coordination, strength and range of motion. Using all four extremities helps with weakness in any area, even general weakness, and can positively affect tone.

Foot trykes were developed in response to requests from therapists for a traditional tricycle for riders with special needs.

Hand trykes are designed for persons whose lower limbs lack function or those who need special therapy for the upper extremity.

Step 3: Take rider's measurements from the front of this Assessment Form and compare them to our Sizing Chart. (This will narrow the choices considerably)

Step 4: Choose a drive. (If it doesn't mention a choice, then ignore this step.)

Amtrykes come with two drive possibilities: **fixed** drive or **geared** drive. Tryke models have been carefully designed so the drive matches the therapeutic goals of the equipment; therefore all drives are not available on all trykes.

A **fixed** drive, commonly known as a 'fixie,' works on a mechanical level to help individuals make a full pedal rotation. The foot crank is constantly in motion for full therapeutic effect. Coasting is not possible with a positive drive; when limb motion stops, the bike does as well.

Riders of **geared** drive trykes must be able to make a complete pedal rotation on their own. This tryke is suitable for riders who need help with balance and a stable sitting position and who have the cognitive ability to steer successfully and apply the coaster brake or hand brake. The key feature of geared trykes is the ability to coast.

Step 5: Chose any adaptations and/or accessories needed by the rider. *Each tryke model can only be customized in the ways noted in its own model section in the Tryke Selection Forms or with the generic accessories listed below. If a customization option is not listed, it is because of design or other practical constraints in offering it on a particular model.*

GENERIC ACCESSORIES (not model specific)

Fun Items	<input type="checkbox"/> License Plate	<input type="checkbox"/> Water Bottle w Cage	<input type="checkbox"/> Mirror
Leg & Foot Items	<input type="checkbox"/> Foot Cups (pair): <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Knee Adductor Strap: <input type="checkbox"/> Small <input type="checkbox"/> Large <input type="checkbox"/> Medium	<input type="checkbox"/> Pedal Block (1 = 3/4") _____ qty
Hand Items	<input type="checkbox"/> Variable Range of Motion Kit <i>(only for Hand & Foot Cycles)</i>	<input type="checkbox"/> Wrist Wraps <i>(Includes right & left)</i> <input type="checkbox"/> X-Small <input type="checkbox"/> Large <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Wrist Brace Mitt: <input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> XX-Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Small <input type="checkbox"/> Small

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name: _____

FOOT TRYKES

1410 FOOT TRYKE – (10" wheels, fixed drive) 50-FC-0100 Arm Length 12-20", Leg Length 15-21"

NEW: Standard with Rear Steering (can convert to push bar) and safety flag

Option #1 – Standard Seating System: Blue Bucket Seat

Option #2 – Alternate Seating System: Snappy Seat System

Snappy Seat Accessories: Notched Seat Bottom Cushion* Laterals Head Rest
 Notched Wide Bottom Cushion* XL Seat Back**

Option #3 – Alternate Seating System: OLD Blue Bucket Seat**

(Less expensive option for very small riders. Cannot be used with an H-Harness or Separator Cube.)

1410 Accessories: Separator Cube Pull/Steering Bar H-Harness - 11.5"
 Heavy-Duty Push Bar – for use with notched seats

1412 ProSeries FOOT TRYKE (12" wheels, fixed drive) – 50-FC-1412 Arm Length 16-24", Leg Length 23-28"

NEW: Standard with Rear Steering (can convert to push bar), Medium Pommel Saddle Seat and safety flag

Option #1 – Standard Seating System: Medium Pommel Saddle Seat/1600 Simple Seatback

Seat Bottom Alternates: Bench Seat Large Pommel Saddle Seat Saddle Seat Skinny Saddle Seat (8.5x6")

Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System (push grip, backpad, med. pommel & 2 laterals)

Seat Bottom Alternates: Bench Seat Large Pommel Saddle Seat Saddle Seat Skinny Saddle Seat (8.5x6")
Trunk Support Options: ProSeries Full Padded Back ProSeries Head Rest ProSeries Lumbar Pad
 Recumbent 10 Degree Seat Post

Option #3 – Alternate Seating System: Gray Bucket Seat

Option #4 – Alternate Seating System: Snappy Seat System *(too small for many 1412 riders; check seat measurements.)*

Snappy Seat Accessories: Notched Seat Bottom Cushion* Laterals Head Rest
 Notched Wide Bottom Cushion* XL Seat Back** Separator Cube

1412 Accessories: Pull/Steering Bar Pedal Toe Pulley H-Harness - 11.5"
 Solid Tires Front Hand Brake Kit (cannot be used with Rear Steering)
 ½" Expanding Pedals Hi-Rise Handlebars: 5" (BMX) 8" ½" Exercise Pedals
 Heavy-Duty Push Bar Knee Separator: 2" 5" 9" Extender Tube

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

NEW! High Functioning Set Up (Both Options Allow Coasting) Can be used with 1416, 1420, 1420XL

(Rider must have the ability to steer, apply brake and make a complete pedal rotation on their own.)

Comes Standard with Pommel Saddle Seat (1412 – M; 1416 & 1420 – L), 1600 Simple Seat Back, Hi-Rise Handlebars, Exercise pedals and safety flag. **MUST** choose Option 1 or 2.

Choose Tryke: 1416 1420 1420XL

Option #1: (switch to existing freewheel hub and add Disk Brake Kit - necessary for safety) Adds the ability to coast

Dual Hand Brake (Option 1 only)

Option #2: (3-Speed Kit and coaster brake) Adds the ability to coast and switch gears for hilly terrain – Replaces the 1600 Series

Standard Seating System (for either option): Pommel Saddle Seat/1600 Simple Seat Back

Seat Bottom Alternates: Bench Seat Tractor Seat with bracket Regular Saddle (11x9")
 Medium Saddle (10.5x10") Large Saddle (13x12")

Adaptive Accessories: ½" Expanding Pedals Telescoping Loop Handlebars 1400 Seat Back Set (w/ laterals)

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name:

FOOT TRYKES

1416 ProSeries FOOT TRYKE (16" wheels, fixed or freewheel drive) – 50-FC-1416 Arm Length 16-24", Leg Length 26-32"

NEW: Standard with Medium Pommel Saddle Seat and safety flag

Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, med. pommel & 2 laterals)

Seat Bottom Alternates: Bench Seat Tractor Seat with bracket Large Pommel Saddle Seat Saddle Seat

Trunk Support Options: ProSeries Full Padded Back ProSeries Head Rest ProSeries Lumbar Pad
 Recumbent 10 Degree Seat Post

Option #2 – Alternate Seating System: Medium Pommel Saddle Seat/1600 Simple Seat Back

Seat Bottom Alternates: Bench Seat Tractor Seat with bracket Large Pommel Saddle Seat Saddle Seat

Option #3 – Alternate Seating System: Gray Bucket Seat

1416 Accessories: Pull/Steering Bar Pedal Toe Pulley H-Harness - 11.5"
 Rear Steering Kit Calf & Leg Supports ½" Exercise Pedals
 Hi-Rise Handlebars: 8" 12" Solid Tires
 Heavy-Duty Push Bar Knee Separator: 2" 5" 9" extender tube

1420 ProSeries FOOT TRYKE (20" wheels, fixed or freewheel drive) – 50-FC-1420 Arm Length 20-30", Leg Length 29-35"

NEW: Standard with Large Pommel Saddle Seat and safety flag

Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, large pommel & 2 laterals)

Seat Bottom Alternates: Bench Seat Tractor Seat w/ bracket Medium Pommel Saddle Seat Saddle Seat

Trunk Support Options: ProSeries Head Rest ProSeries Full Padded Back ProSeries Lumbar Pad
 Recumbent 10 Degree Seat Post

Option #2 – Alternate Seating System: Large Pommel Saddle Seat/1600 Simple Seatback

Seat Bottom Alternates: Bench Seat Tractor Seat w/ bracket Medium Pommel Saddle Seat Saddle Seat

Option #3 – Alternate Seating System: Gray Bucket Seat

1420 Accessories: Pull/Steering Bar Pedal Toe Pulley H-Harness - 11.5"
 Rear Steering Kit Calf & Leg Supports ½" Exercise Pedals
 12" Hi-Rise Handlebars Solid Tires
 Heavy-Duty Push Bar Knee Separator: 2" 5" 9" extender tube

1420XL ProSeries FOOT TRYKE (20" wheels, fixed/freewheel drive) – 50-FC-1420-XL Arm Lgth 22-34", Leg Lgth 30-45"

NEW: Standard with Large Pommel Saddle Seat and safety flag

Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, large pommel & 2 laterals)

Seat Bottom Alternates: Bench Seat Tractor Seat w/ bracket Medium Pommel Saddle Seat Saddle Seat

Trunk Support Options: ProSeries Head Rest ProSeries Full Padded Back ProSeries Lumbar Pad
 Recumbent 10 Degree Seat Post

Option #2 – Alternate Seating System: Large Pommel Saddle Seat/1600 Simple Seat Back

Seat Bottom Alternates: Bench Seat Tractor Seat w/ bracket Medium Pommel Saddle Seat Saddle Seat

Option #3 – Alternate Seating System: Large Wheelchair Seat

Wheelchair Seat Accessories: Wheelchair Seat Bar Ends Swing Away Arms

1420 XL Accessories: Pull/Steering Bar Toe Pedal Pulley H-Harness - 11.5"
 Rear Steering Kit Calf & Leg Supports ½" Exercise Pedals
 Hi-Rise Handlebars 12" 23"
 Heavy-Duty Push Bar Knee Separator: 2" 5" 9" extender tube

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name: _____

FOOT TRYKES

<input type="checkbox"/> JT-2000 FOOT TRYKE (Recumbent 14-Speed) – 50-FC-2000 (advanced riders only) Arm Length 20-28", Leg Length 30-41"				
JT-2000 Accessories:	<input type="checkbox"/> Exercise Pedals	<input type="checkbox"/> XL Exercise Pedals	<input type="checkbox"/> Expanding Pedals	<input type="checkbox"/> Digital Speedometer
	<input type="checkbox"/> Dual Hand Brake	<input type="checkbox"/> Toe Clips	<input type="checkbox"/> Basket	<input type="checkbox"/> Rearview Mirror

<input type="checkbox"/> JT-2300-USS FOOT TRYKE (Recumbent 14-Speed with Under Seat Steering) – 50-FC-2300-USS (advanced riders only) Arm Length 20-28", Leg Length 30-41"				
JT-2300-USS Accessories:	<input type="checkbox"/> Exercise Pedals	<input type="checkbox"/> XL Exercise Pedals	<input type="checkbox"/> Expanding Pedals	<input type="checkbox"/> Digital Speedometer
	<input type="checkbox"/> Dual Hand Brake	<input type="checkbox"/> Toe Clips	<input type="checkbox"/> Basket	<input type="checkbox"/> Rearview Mirror

<input type="checkbox"/> TP-3000 TADPOLE FOOT TRYKE (Recumbent 21-Speed) – 50-FC-3000 (advanced riders only) Arm Length 17-29", Leg Length 36-45"				
Tadpole Accessories:	<input type="checkbox"/> Exercise Pedals	<input type="checkbox"/> XL Exercise Pedals	<input type="checkbox"/> Expanding Pedals	
	<input type="checkbox"/> Dual Hand Brake	<input type="checkbox"/> Toe Clips	<input type="checkbox"/> Rearview Mirror	

Therapist Assembly Notes or Comments (if any):

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name: _____

HAND TRYKES

<input type="checkbox"/> AM-10 HAND TRYKE (12" wheels, 2.5" crank arms & permanent foot platform) 50-HFC-0105 Arm length 13-17" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**
AM-10 Accessories: <input type="checkbox"/> Vertical Hand Grips <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Heavy Duty Push Bar (<i>notched seats</i>) <input type="checkbox"/> Separator Cube

<input type="checkbox"/> AM-12Small HAND TRYKE (12" wheels, 3" crank arms & permanent foot platform) 50-HFC-0110 Arm length 15-20" New: Now standard with Saddle Seat, Rear Steering (converts to push bar by inserting pin) and safety flag.
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat & 1600 Simple Seat Back
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**
<input type="checkbox"/> Option #3 – Alternate Seating System: Blue Bucket Seat
AM-12S Accessories: <input type="checkbox"/> Vertical Hand Grips <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Heavy Duty Push Bar (<i>notched seats</i>) <input type="checkbox"/> Separator Cube

<input type="checkbox"/> AM-12 HAND TRYKE – (12" wheels, 4" crank arms & permanent foot platform) 50-HC-0210 Arm Length 14-23" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat
<input type="checkbox"/> Option #4 – Alternate Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Separator Cube
AM-12 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube <input type="checkbox"/> Vertical Hand Grips

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

<input type="checkbox"/> AM-16 HAND TRYKE (16" wheels, 5" crank arms & permanent foot platform) 50-HFC-0411 Arm length 18-27" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Padded Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat
AM-16 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube <input type="checkbox"/> Vertical Hand Grips

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name:

HAND TRYKES

NEW: 1020 "Junior" HAND TRYKE (20" wheels & 3-Speed drive train) 50-HC-1020 Arm length 19-30"

Option #1 – Standard Seating: Small Wheelchair Seat (3.5" narrower in width and depth than the Large)

Seat Alternates: Large Wheelchair Seat

1020 Accessories: Wheelchair Seat Bar Ends Swing-Away Arms Rearview Mirror

1024 HAND TRYKE – (24" wheels & 3-Speed drive train) 50-HC-1024 Arm length 22-26"

Option #1 – Standard Seating: Large Wheelchair Seat (3.5" broader in width and depth than the Small)

Seat Alternates: Small Wheelchair Seat

1024 Accessories: Wheelchair Seat Bar Ends Swing-Away Arms Rearview Mirror

HP-1000 RECREATION HAND TRYKE – (16" front wheel & 7-Speed grip shift) 50-HC-1000 Arm length 22-26"

HP-1000 Accessories: Wheelchair Seat Bar Ends Swing-Away Arms Rearview Mirror

Therapist Assembly Notes or Comments (if any):

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name:

HAND-FOOT TRYKES

AM-10 HAND-FOOT TRYKE (10" wheels & 2.5" crank arms) 50-HFC-0105 Arm length 13-17", leg length 15-21"

New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.

Option #1 – Standard Seating System: Blue Bucket Seat

Option #2 – Upgraded Seating System: Snappy Seat System

Snappy Seat Accessories: Notched Seat Bottom Cushion* Laterals Head Rest

Notched Wide Bottom Cushion* XL Seat Back**

AM-10 Accessories: AM Toe Pulley Pull/Steering Bar H-Harness - 11.5"

Heavy-Duty Push Bar (notched seats) Vertical Hand Grips Separator Cube

AM-12Small HAND-FOOT TRYKE (12" wheels & 3" crank arms) 50-HFC-0110 Arm length 15-20", leg length 19-24"

New: Now standard with Saddle Seat, Rear Steering (converts to push bar by inserting pin) and safety flag.

Option #1 – Standard Seating System: Saddle Seat & 1600 Simple Seat Back

Option #2 – Upgraded Seating System: Snappy Seat System

Snappy Seat Accessories: Notched Seat Bottom Cushion* Laterals Head Rest

Notched Wide Bottom Cushion* XL Seat Back**

Option #3 – Upgraded Seating System: Blue Bucket Seat

AM-12S Accessories: AM Pedal Toe Pulley Pull/Steering Bar H-Harness - 11.5"

Heavy-Duty Push Bar (notched seats) Vertical Hand Grips Separator Cube

9/16" Exercise pedals Solid Rear Tires

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

AM-12 HAND-FOOT TRYKE (12" wheels & 4" crank arms) 50-HFC-0210 Arm length 14-23", leg length 21-29"

New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.

Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back

Seat Bottom Alternates: Bench Seat Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Option #2 – Upgraded Seating System: ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)

Seat Bottom Alternates: Bench Seat Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Trunk Support Upgrades: ProSeries Full Padded Back ProSeries Head Rest ProSeries Lumbar Pad

Option #3 – Upgraded Seating System: Gray Bucket Seat

Option #4 – Upgraded Seating System: Snappy Seat System (too small for many AM-12 riders; check seat measurements.)

Snappy Seat Accessories: Notched Seat Bottom Cushion* Laterals Head Rest

Notched Wide Bottom Cushion* XL Seat Back** Separator Cube

AM-12 Accessories: AM Pedal Toe Pulley Pull/Steering Bar H-Harness - 11.5"

Knee Separator: 2" 5" 9" Extender Tube Vertical Hand Grips

AM Leg Rests 9/16" Exercise pedals Solid Rear Tires

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR

Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Recipient's Name:

AM-16 HAND-FOOT TRYKE (16" wheels & 5" crank arms) 50-HFC-0411 Arm length 18-27", leg length 24-36"

New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.

Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back

Seat Bottom Alternates: Bench Seat Tractor Seat Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)

Seat Bottom Upgrades: Bench Seat Tractor Seat Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Trunk Support Upgrades: ProSeries Full Padded Back ProSeries Head Rest ProSeries Lumbar Pad

Option #3 – Upgraded Seating System: Gray Bucket Seat

AM-16 Accessories: AM Pedal Toe Pulley Pull/Steering Bar H-Harness - 11.5"

Knee Separator: 2" 5" 9" Extender Tube AM Leg Rests XL Exercise pedals

Solid Rear Tires 9/16" Exercise pedals Vertical Hand Grips

~~**AM-20 HAND-FOOT TRYKE**~~ (Recumbent Single Speed) – 50-HFC-0610 Arm length 18-30", leg length 32-46"

Tryke is currently being used **ONLY** for research.

AM-20 Accessories: Exercise Pedals XL Exercise Pedals Expanding Pedals

Wheelchair Seat Bar Ends Swing Away Arms Calf & Leg Supports (Expanding pedals)

Therapist Assembly Notes or Comments (if any):

Please scan & email application & photo to:
applications@capitalcityambucs.org

OR


Mail to: Capital City AMBUCS
3207 Dancy St. Austin, Tx 78722

Amtryke Sizing Chart

TRYKE TYPE (How will the tryke be propelled?)	RIDER LEG LENGTH (Inches from center of hip to bottom of shoe.)	RIDER ARM LENGTH (Inches from middle of shoulder to center of digit crease.)	MODEL	RIDER WEIGHT	RIDER MAX HEIGHT	TRYKE WEIGHT	WHEEL SIZE	TTRYKE HEIGHT	TRYKE LEGNTH	TRYKE WIDTH
				(Pounds)	(Inches)	(Pounds)	(Inches)	(Inches)	(Inches)	
Hand & Foot	15-21	13-17	AM-10	55	40	45	10	24	38	21
	19-24	15-20	AM-12S	150	40	45	12	27	38	24
	21-29	14-23	AM-12	150	47	45	12	36	60	32
	24-36	18-27	AM-16	175	66	55	16	36	68	33
	32-46	18-30	AM-20	FOR RESEARCH ONLY						
Foot	15-21	12-20	1410	55	40	45	10	24	38	21
	23-28	16-24	1412	125	42	72	12	13	43	27
	26-32	16-24	1416	175	60	74	16	49	58	30
	29-35	20-30	1420	250	68	74	20	50	64	30
	30-45	22-34	1420XL	275	76	89	20	43	72	29
	30-41	20-28	JT-2000/ JT2300USS	250	74	80	20	48	72	32
	36-45	17-29	TP-3000	300	75	47	20	33	63-70	33
	up to 37	19-30	1020	250	67	74	20	41	69	30
	up to 41	22-26	1024	250	72	85	24	45	75	32
	up to 41	22-26	HP-1000	250	74	83	Front 16 Rear 20	45	72	32

All trykes in the Hand & Foot section can be converted to Hand trykes – except the AM-20.

- A Center of Shoulder
- B Center of Elbow
- C Center of Digit Crease
- D Center of Hip (greater trochanter)
- E Center of Knee
- F Bottom of Foot



RIDER'S MEASUREMENTS	
Arm Measurements (Inches) Total Length	
Left	A to B + B to C = _____
Right	A to B + B to C = _____
Trunk	A to D = _____
Leg Measurements (Inches) Total Length	
Left	D to E + E to F = _____
Right	D to E + E to F = _____

HELMET SIZING	
Youth Sizes	Head Circumference Inches
XXS	18.5 to 19.5
XS	20.5 to 22
S/M	22 to 23.6
L/XL	23.6 to 25.75
Adult Sizes	Head Circumference Inches
S/M	22 to 23.6
L/XL	23.6 to 25.75